


- ▲ 8 forceps in a custom titanium tray
- ▲ Designed by Dr. Robert Spetzler for the unique challenges of microsurgery
- ▲ Crisp edges and polished surfaces
- ▲ Color coded handles keep instruments organized in the custom tray



N23105-SET




#1 Straight Forceps
0.5mm Tips, 21cm
N4401




#5 Curved Up Forceps
0.3mm Tip with Interlocking
1x2 Teeth, 21cm
N4405



#2 Straight Forceps
0.5mm Tip with Interlocking
1x2 Teeth, 18cm
N4402



#6 Curved Down Forceps
0.3mm Tip with Interlocking
1x2 Teeth, 21cm
N4406



#3 Straight Forceps
0.5mm Tip with Interlocking
1x2 Teeth, 21cm
N4403



#7 Ring Forceps
2.0mm Ring Tip, 21cm
N4407



#4 Straight Forceps
0.3mm Tip with Interlocking
1x2 Teeth, 21cm
N4404



#8 Cup Forceps
1.5mm Cup Tip, 21cm
N4408

The leader in surgical instrument innovation...